

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody (orb1084062)

Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein (1-419/419 aa)
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.05% Proclin300.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.