

Product Datasheet

Sulf2 peptide (orb108407)

Description	Sulf2 blocking peptide
Conjugation	Unconjugated
Tested Applications	FA
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Blocking peptide for use with Sulf-2 antibodies (Cat. Nos. X1858P & X1865P).
Uniprot ID	Q8IWU5
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.