

Product Datasheet

SIAH1 peptide (orb108408)

Description	SIAH1, E3 Ubiquitin Ligase blocking peptide
Conjugation	Unconjugated
Tested Applications	FA
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Blocking peptide for use with SIAH1 antibodies (Cat. Nos. X1859P & X1864P).
Uniprot ID	Q8IUQ4
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.