

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody (orb1084086)

Description

SARS-CoV-2 (2019-nCoV)

Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from SARS-CoV-2 Nucleocapsid Protein
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.