

Product Datasheet

SARS-CoV-2 (2019-nCoV) Spike protein S2 antibody (orb1084091)

Description	SARS-CoV-2 (2019-nCoV) Spike protein S2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from SARS-CoV-2 Spike Protein (1091-1190/1273 aa)
Target	SARS-CoV-2 (2019-nCoV) Spike protein S2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	141/65 kDa
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.