

Product Datasheet

Lass6 peptide (orb108433)

Description	Lass6 blocking peptide
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide Immunogen
Preservatives	Provided as solution in phosphate buffered saline
Form/Appearance	Provided as solution in phosphate buffered saline
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	For use with Lass6 polyclonal antibodies (Cat. No. X2301P & X2303P).
Uniprot ID	Q6ZMG9
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.