

Product Datasheet

Roundabout homolog 1 peptide (orb108435)

Description	Roundabout homolog 1, ROBO1 blocking peptide
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide derived from human ROBO1 protein
Preservatives	Provided as solution in phosphate buffered saline
Form/Appearance	Provided as solution in phosphate buffered saline
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	For use with ROBO1 polyclonal antibodies (Cat. No. X2313P and X2334P).
Uniprot ID	Q9Y6N7
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.