



Product Datasheet

HSP60 antibody (orb10846)


Description	HSP60 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human HSP60 (501-573/573 aa)
Target	HSP60
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P10809
Dilution Range	WB=1:500-2000, ICC=1:100, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of HepG2 cell using HSP60 p...



Western blot:Kidney (Mouse);Primary:Anti...



Western blot:HeLa(Human) Cell;Primary:An...