

Product Datasheet

PINK1 peptide (orb108486)

Description	PINK1 (E417G Mutant) blocking peptide
Conjugation	Unconjugated
Tested Applications	Blocking
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Optimal concentration should be evaluated by serial dilutions.
Uniprot ID	Q9BXM7
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.