

Product Datasheet

Bovine Collagen I protein (orb108488)

Description	Bovine Collagen I protein
Species/Host	Bovine
Reactivity	Bovine
Conjugation	FITC
Immunogen	Bovine skin
Preservatives	Provided as a 1 mg/ml solution in 0.01M acetic acid.
Form/Appearance	Provided as a 1 mg/ml solution in 0.01M acetic acid.
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Used for the assay of mammalian collagenase activity. To heat denature collagen for preparing a substrate for gelatinase, transfer collagen solution to a brown glass bottle and heat to 80°C.
Clonality	Monoclonal
Purity	> 95%
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	6 months from date of receipt.