

Product Datasheet

Anti-Phospho-MEK7 (Ser271 + Thr275)/Cy7 (orb1084928)

Description	Anti-Phospho-MEK7 (Ser271 + Thr275)/Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy7
Tested Applications	FC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MKK7 around the phosphorylation site of Ser271/Thr275
Target	Phospho-MEK7 (Ser271 + Thr275)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14733
Dilution Range	Flow-Cyt=1µg /Test,IF=1:100-500
Expiration Date	12 months from date of receipt.