


Product Datasheet

E1A antibody (orb108521)

Description	Sheep polyclonal to E1A
Species/Host	Sheep
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Recombinant E1A protein recognizing a 43kD protein
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Western Blot suggested dilution 1- 10 ug /ml
Isotype	IgG
Uniprot ID	P03254
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunoprecipitation
analysis of HCT116
C...