

## Product Datasheet

### BCL2 antibody (orb108524)

<b>Description</b>	Mouse monoclonal antibody to BCL2
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthetic peptide corresponding to amino acids 3 to 14 of the human bcl-xl sequence
<b>Preservatives</b>	Provided as solution in phosphate buffered saline with no preservative or carrier protein added.
<b>Form/Appearance</b>	Provided as solution in phosphate buffered saline with no preservative or carrier protein added.
<b>Concentration</b>	See vial for concentration
<b>Storage</b>	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Application: Western Blot at 1-5 ug/ml.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>Q07817</b>
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Cell lysates us...