

Product Datasheet

Caspase 3 antibody (orb108527)

Description	Rabbit polyclonal antibody to Caspase 3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length recombinant caspase-3 protein
Preservatives	Provided as solution in phosphate buffered saline, pH 7.3, with 1.0 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline, pH 7.3, with 1.0 mg/ml rAlbumin and 0.05% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: western blot, use at 2-10 µg/ml
Isotype	IgG
Uniprot ID	P42574
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.