

## Product Datasheet

### BCL2 antibody (orb108529)

<b>Description</b>	Provided as solution in phosphate buffered saline with 0.08% sodium azide
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Concentration</b>	See vial for concentration
<b>Storage</b>	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Western Blots at 1-10 ug/ml
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mice were immunized with a synthetic peptide corresponding to amino acids 61 to 76 of the mouse bcl-2 sequence
<b>Uniprot ID</b>	<b>P10415, P49950</b>
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
<b>Expiration Date</b>	12 months from date of receipt.