

Product Datasheet

E1A antibody (orb108532)

Description	Mouse monoclonal to E1A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Immunoprecipitation (1-2 ug/ml) and 2-5 ug/ml for Immunofluorescence.
Isotype	IgG2a
Clonality	Monoclonal
Source	The immunogen for the anti-2 E1A was full length recombinant adenovirus 2 E1A protein.
Uniprot ID	P03254
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.