

Product Datasheet

Cyclin Dependant Kinase 7 antibody (orb108544)

Description	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: 1-10 ug/ml for Western Blot, Immunohistochemistry, Immunoprecipitation
Isotype	IgG2b
Clonality	Monoclonal
Source	Rabbits were immunized with a recombinant 221 amino acid fragment of the C terminus of p40mo15 recognizing a 40kD protein
Uniprot ID	P50613
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.