


Product Datasheet

EBNA 3A antibody (orb108558)

Description	Sheep polyclonal to EBNA 3A
Species/Host	Sheep
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Western blotting at 0.25-1.00 µg/ml, Dot blotting, ELISA, Immunocytochemistry and Immunohistochemistry on frozen and paraffin embedded tissues. Optimal concentrations must be determined by titration.
Isotype	IgG
Source	Sheep were immunized with recombinant EBNA-3A protein
Uniprot ID	Q69138
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of cell lines using...