

Product Datasheet

MAGI antibody (orb108564)

Description	Sheep polyclonal to MAGI
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein corresponding to PDZ domains 4-5 of MAGI-3.
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody detects MAGI-3 by Western blot at 1-10 µg/ml. Detects a 130 kDa band in cells transfected with full-length MAGI-3 protein. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Uniprot ID	Q9H2V6
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of cell lysates us...