

Product Datasheet

p38 HOG Kinase antibody (orb108568)

Description	Sheep polyclonal to p38 HOG Kinase
Species/Host	Sheep
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Full length p38 HOG kinase protein
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: 1-10 ug/ml for immunoprecipitation of the catalytically active form of p38HOG kinase in vitro kinase assay.
Isotype	IgG
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.