

Product Datasheet

TGF beta II receptor antibody (orb108571)

Description	Sheep polyclonal to TGF beta II receptor
Species/Host	Sheep
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A recombinant protein containing chick transforming growth factor type II receptor a.a.27-85
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: Western Blot (1:100), Immunohistology, Empirical
Isotype	IgG
Uniprot ID	Q90999
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.