

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleoprotein ELISA Kit (orb1085810)

Description	SARS-CoV-2 (2019-nCoV) Nucleoprotein ELISA Kit
Range	0.78 - 25ng/ml
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates and cell culture extracts
Sensitivity	0.4 ng/ml
Expiration Date	Please enquire.