

Product Datasheet

p130 antibody (orb108605)

Description	Monoclonal antibody provided as solution in phosphate buffered saline with 0.08% sodium azide
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Western Blots at 1:1000
Isotype	IgG1
Clonality	Monoclonal
Source	The immunogen used was a recombinant protein to the N terminus of p130.
Uniprot ID	Q64700 , O55081 , Q08999
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.