


Product Datasheet

c-myc antibody (orb108607)

Description	Mouse monoclonal to c-myc
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: Western Blot, Immunoprecipitation (Both with and without SDS) and immunohistochemistry (see image)
Isotype	IgG2a
Clonality	Monoclonal
Source	Hybridoma produced by the fusion of splenocytes from mice immunized with recombinant full-length human c-Myc protein and mouse myeloma cells.
Uniprot ID	P01106, P01108, P09416
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of normal k...