

Product Datasheet

Human CCL24/Eotaxin-2 ELISA KIT (orb1086171)

Description	Human CCL24/Eotaxin-2 ELISA KIT
Reactivity	Human
Range	31.20 - 2,000 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	15 pg/mL
Expiration Date	Please enquire.