

## Product Datasheet

### Human CD27/TNFRSF7 ELISA KIT (orb1086189)

|                        |                                                                                                              |
|------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Human CD27/TNFRSF7 ELISA KIT                                                                                 |
| <b>Reactivity</b>      | Human                                                                                                        |
| <b>Range</b>           | 125-8000 pg/mL                                                                                               |
| <b>Storage</b>         | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.) |
| <b>Note</b>            | For research use only                                                                                        |
| <b>Sample Types</b>    | cell culture supernates, serum, and plasma.                                                                  |
| <b>Sensitivity</b>     | 62 pg/mL                                                                                                     |
| <b>Expiration Date</b> | Please enquire.                                                                                              |