



Product Datasheet

IASPP antibody (orb10862)


Description	IASPP antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human IASPP (15-100/828 aa)
Target	IASPP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89 kDa
Uniprot ID	Q8WUF5
Dilution Range	ICC=1:100, Flow-Cyt=1µg/Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P of keloid tissue
(IASPP antibody a...)



Western blot analysis
of IASPP expressio...



Immunohistochemical
analysis of paraffin...