

Product Datasheet

Protein Phosphatase 1 gamma 1 antibody (orb108620)

Description

Rabbit polyclonal to Protein Phosphatase 1 gamma

Species/Host

Rabbit

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IP, WB

Applications

Western blot analysis of rat brain homog...

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Western blotting (5-10 ug/ml) and immunoprecipitation.

Isotype

IgG

Source

The immunogen used to generate the purified antibody was a peptide conjugated to KLH corresponding to the sequence NH₂-Thr-Pro-Pro-Arg-Gly-Met-Ile-Thr-Lys-Gln-Ala-Lys-Lys-COOH. This peptide antibody corresponds to C-terminal peptide (amino acids 311-323) of PP1 gamma 1 catalytic subunit having a MW of 37kDa.

Uniprot ID

P37139, **P36873**

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.