

Product Datasheet

Mouse Complement Component C5a ELISA KIT (orb1086277)

Description	Mouse Complement Component C5a ELISA KIT
Reactivity	Mouse
Range	15.60-1000 pmol/l
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	7 pmol/l
Expiration Date	Please enquire.