

Product Datasheet

Protein Phosphatase 2 A/B delta antibody (orb108628)

Description

Rabbit polyclonal to Protein

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide conjugated to KLH derived from the N-terminus of PP2A/B delta
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Detects PP2A/B delta protein in rat testis by Western blot. Detects a 53 kDa band in rat brain, heart, lung, placenta, skeletal muscle, spleen, testis and thymus tissues lysates. Optimal concentration should be evaluated by serial dilutions. Applications: Western Blotting, Immunofluorescence
Isotype	IgG
Uniprot ID	Q925E7 , P56932 , Q14738
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.