

Product Datasheet

Mouse CCL11/Eotaxin-1 ELISA KIT (orb1086283)

Description	Mouse CCL11/Eotaxin-1 ELISA KIT
Reactivity	Mouse
Range	7.81 - 500 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	4 pg/mL
Expiration Date	Please enquire.