


Product Datasheet

Protein Phosphatase 2 A/B antibody (orb108629)

Description	This antibody is specific to pp2A/B isoforms.
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Application: Immunoblotting 1:500 to 1:1,000.
Isotype	IgG
Source	The immunogen used to generate the purified antibody was a peptide conjugated to KLH corresponding to the sequence NH ₂ -Val-Ser-Ser-Pro-His-Phe-Gln-Val-Ala-GluArg-Ala-Leu-Tyr-COOH.
Uniprot ID	N/A
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain extra...