


Product Datasheet

Vav-2 SHD antibody (orb108647)

Description	Sheep polyclonal to Vav-2 SHD
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A his tagged fusion protein Vav-2 corresponding to amino acid sequence 578 - 878 and having a molecular weight of 100kD
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: Western Blot 1:5,000. Immunoprecipitation -Empirical
Isotype	IgG
Uniprot ID	P52735
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of HEK 283

...