


Product Datasheet

LEF/TCF/HMG antibody (orb108656)

Description	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for Western blotting (1-5 µg/ml) and immunofluorescence on paraformaldehyde-fixed cells (1:1000 or 1 µg/ml). Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG1
Clonality	Monoclonal
Source	The immunogen for the anti-LEF-1 antibody was a full length protein.
Uniprot ID	Q9UJU2
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Cos-1 cel...