

Product Datasheet

Calpain I, II antibody (orb108667)

Description	Provided as a concentration of 1 mg/ml in 0.02M sodium phosphate, pH 7.5, 0.15M sodium chloride, 50%...
Species/Host	Mouse
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Use at 0.5-1 ug/ml for ELISA and western blot, 1-2ug/ml for immunoprecipitation.
Isotype	IgG1
Clonality	Monoclonal
Source	Hybridoma produced by the fusion of splenocytes from mice immunized with human calpain with mouse myeloma.
Uniprot ID	P17655
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.