

Product Datasheet

Phosphotyrosine antibody (FITC) (orb108675)

Description	Provided as sterile solution in 20 mM sodium phosphate, 150 mM sodium chloride, 50% glycerol, 3 mM s...
Species/Host	Mouse
Reactivity	All
Conjugation	FITC
Tested Applications	ELISA, IHC-P, IP, WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	A dilution of 5-10 ug/ml is suggested for immunohistochemistry or immunofluorescence. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG2b
Clonality	Monoclonal
Source	Hybridoma produced from Balb/C mice immunized with phosphotyrosine coupled to carrier protein.
Uniprot ID	N/A
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.