

Product Datasheet

phosphotyrosine antibody (Peroxidase) (orb108677)

Description

Mouse monoclonal antibody to

Peroxidase.

Species/Host

Mouse

Reactivity

All

Conjugation

HRP

Tested

ELISA, IHC-P, IP, WB

Applications
Immunogen

Hybridoma produced from Balb/C mice immunized with phosphotyrosine coupled to carrier protein.

Preservatives

Provided in 60% glycerol, phosphate buffered saline, and 0.1% 2-Chloroacetamide

Form/Appearance

Provided in 60% glycerol, phosphate buffered saline, and 0.1% 2-Chloroacetamide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Application: A dilution of 1:2,000 is suggested for western blot or ELISA, 1:500 for immunoprecipitation, immunohistochemistry, or immunocytochemistry.

Isotype

IgG2b

Clonality

Monoclonal

Uniprot ID

N/A

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.