

Product Datasheet

phosphotyrosine antibody (Peroxidase) (orb108679)

Description	Provided in 60% glycerol, phosphate buffered saline, and 0.1% 2-Chloroacetamide
Species/Host	Mouse
Reactivity	All
Conjugation	HRP
Tested Applications	ELISA, IHC-P, IP, WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Working dilution of 1:2,000 is suggested for western blot or ELISA, 1:500 for immunohistochemistry or immunocytochemistry.
Isotype	IgG2b
Clonality	Monoclonal
Source	Hybridoma produced from Balb/C mice immunized with phosphotyrosine coupled to carrier protein.
Uniprot ID	N/A
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.