


Product Datasheet

ERBB4 antibody (orb10868)

Description	Rabbit polyclonal antibody to ERBB4
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human ErbB4/HER4(NP_005226.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	147kDa
Uniprot ID	Q15303
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of rat...