


Product Datasheet

BrdU antibody (orb108686)

Description	Mouse monoclonal antibody to BrdU
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P
Immunogen	Bromodeoxyuridine conjugated to Helix Pomatia Haemocyanin used to immunize BALB/c mice
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Reacts with cells containing incorporated BrdU, showing a clear, nucleus confined, speckled pattern. Optimal concentration should be evaluated by serial dilutions. Applications: Flow Cytometry Staining, Immunohistochemistry.
Isotype	IgG1
Clonality	Monoclonal
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse in...