

Product Datasheet

HA agarose antibody (orb108687)

Description	Provided as 2x slurry (0.5 ml of settled gel) containing 20 mM sodium phosphate, 150 mM sodium chlor...
Species/Host	Mouse
Conjugation	Agarose
Tested Applications	IP
Concentration	See vial for concentration
Storage	Product should be stored at 4-8°C. DO NOT FREEZE
Note	For research use only
Application notes	Use 25µl of slurry/100µg of lysate for immunoprecipitation. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG2b
Clonality	Monoclonal
Source	Balb/C mice were immunized with a 9 amino acid petide of influenza virus hemagglutinin
Uniprot ID	Q3TBD2
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.