


Product Datasheet

3BP2 antibody (orb108701)

Description	Sheep polyclonal to 3BP2
Species/Host	Sheep
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Antibody developed using SH2 domain of the 3BP2 protein fused to GST.
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Detects 3BP2 protein by Western blot analysis. Detects a 85 kDa band in murine fibroblasts transfected with 3BP2. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Uniprot ID	Q06649
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant ful...