



Product Datasheet

CD218a antibody (orb1087499)


Description	Mouse monoclonal antibody to CD218a
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, IP, WB
Immunogen	Human CD218a expressed on a transfected NK0 cell line
Target	CD218a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q13478
Entrez	8809
Dilution Range	FC: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
neutrophil
granulocy...



Flow
cytometry
surface
staining
pattern ...