

Product Datasheet

Mouse IFN-alpha / beta R1 Protein (orb1087552)

Description

Mouse IFN-alpha / beta R1, His Tag is expressed

27...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IFN-alpha/beta R1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 47.6 kDa. The protein migrates as 57-66 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse IFN-alpha / beta R1, His Tag (orb1087552) is expressed from human 293 cells (HEK293). It contains AA Glu 27 - Thr 429 (Accession # P33896-1).

MW

47.6 kDa

Uniprot ID


P33896

NCBI

NP_034638.2

Expiration Date

6 months from date of receipt.



Mouse IFN-alpha / beta R1, His Tag (Cat....