


Product Datasheet

Human IL-2 Protein (orb1087555)

Description	Human IL-2 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ala 21 - Thr 1...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 41.9 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-2 Protein, Fc Tag (orb1087555) is expressed from human 293 cells (HEK293). It contains AA Ala 21 - Thr 153 (Accession # P60568-1).
MW	41.9 kDa
Uniprot ID	P60568
NCBI	NP_000577.2
Expiration Date	6 months from date of receipt.



Immobilized
Human IL-2 R
alpha, His
Tag ...