

Product Datasheet

Trypsin Pretreatment Kit, 155 ml (125 ml Buffer 30 ml Enzym) (orb1087604)

Description	Trypsin Pretreatment Kit, 155 ml (125 ml Buffer 30 ml Enzym)
Tested Applications	IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Storage	The solutions should be stored at 2-8°C. Please store the reagents in a dark place and do not freeze them. Under these conditions the solutions are stable up to the expiry date indicated on the label. Do not use product after the expiry date. A positive and a negative control have to be carried out in parallel to the test material. If you observe unusual staining or other deviations from the expected results which could possibly be caused by this reagent, please contact our technical support.
Note	For research use only
Application notes	The components of this Trypsin Pretreatment Kit allow for the preparation of a buffered trypsin solution for enzymatic epitope retrieval.
Expiration Date	Please enquire.