

Product Datasheet

Mouse anti-CD45RO, clone UCHL-1 (monoclonal) (orb1087652)

Description	Mouse anti-CD45RO, clone UCHL-1 (monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tissue culture supernatant with sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	IHC-P (1:200-1:800),
Expiration Date	12 months from date of receipt.