

Product Datasheet

Mouse anti-Carcinoembryonic Antigen (CD66e), clone COL-1 (monoclonal) (orb1087667)

Description	Mouse anti-Carcinoembryonic Antigen (CD66e), clone COL-1 (monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tissue culture supernatant with sodium azide
Concentration	Greater than or equal to 20 mg/L
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	IHC-P (1:100)
Expiration Date	12 months from date of receipt.