


Product Datasheet

chx 10 antibody (orb108770)

Description	Antibody raised against a recombinant protein corresponding to amino acids 264 to 361 derived from t...
Species/Host	Sheep
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Detects Chx10 by Western blot at 0.5 to 1 µg/ml. Detects a 46 kDa band in mouse and rat retinal tissue lysates. Optimal concentration should be evaluated by serial dilutions. Suitable for use against recombinant proteins conjugated to OVA, GST, His tags and other.
Isotype	IgG
Uniprot ID	P58304
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human me...