

Product Datasheet

Mouse anti-Myosin Heavy Chain (slow), clone WB-MHCs (monoclonal) (orb1087715)

Description	Mouse anti-Myosin Heavy Chain (slow), clone WB-MHCs (monoclonal)
Species/Host	Mouse
Reactivity	Canine, Goat, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized tissue culture supernatant containing 15 mM sodium azide.
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-fr (1:20-1:80)
Expiration Date	12 months from date of receipt.